eversense 365
Continuous Glucose Monitoring System


DMS Pro **Quick Set Up**

This setup guide provides a shortcut for establishing an Eversense DMS Pro account tailored for your clinic.
For more instructions, please refer to the DMS Pro User Guide **(https://provider.eversensecgm.com/user-guides/)**.

1 Sign Up for DMS Pro

Go to **https://uspro.eversensedms.com.**
Click on **Sign Up.**

Click on **Sign Up**



---

2 Request an Eversense DMS Pro Account


Complete the form by entering your clinic's information.

- Under **First Name** and **Last Name,** it is recommended to put either your clinic name or the name of the primary healthcare provider associated with the practice. (This is what your patients will see when they share data.)
- Select a **Title** from the dropdown, and complete the remaining required fields.*

Note:
To add individual providers and staff members to your Eversense DMS Pro clinic, see the Eversense DMS Pro user guide.

***You cannot use the same email address for both an Eversense DMS Pro account and a personal Eversense DMS account.**

Click **Submit** when form is completed



---

3 You Will Receive Two Emails

1. Account registration confirmation.
2. Account approval and activation email (may take up to 24 hours). Click on the link† in the email to complete account activation. Then record account information in the box to the right.

†The link expires after 72 hours, but you can request that a new email be sent.

Clinic Name: .....

Email: .....


Password: .....


Clinic ID#: (see Step 4) .....

4 Sign in to the Clinic DMS Pro account

Once logged in, your assigned Clinic ID will appear at the top of home page.


Note:
Patients will need your Clinic ID number when sending a request to share their CGM data with the clinic.






5 Inviting Patients to the Clinic:

- On the left-hand navigation menu, click **Patients,** then select **Invite Patient.**
- Enter the patient's Eversense email (the one they used to create their Eversense account).
- Click **Search Eversense** to locate the patient and send the invitation.








If the patient's email is registered, the Invite patient screen will be displayed. Complete all mandatory information, then hit **Send.** The patient can accept the invitation either via email or in the Eversense App.


Note: If the entered email is not registered with Eversense DMS, you will receive an error message. Confirm with the patient which email address was used when setting up their account.

Until the patient has accepted the invitation, **Pending Patient Consent** will be displayed in the Status column. Once they accept their status will change to **Active**.



6 Accepting Patient Requests

- Patients may request to share their data. Pending requests from patients appear in the left navigation bar. Click on **Patients > Pending Requests.**
- Click **Accept** to add the patient to your clinic. When you accept you will be prompted to assign a provider from your clinic list. Once you have accepted their invitation their status will change to **Active** in your clinic list.





For additional information or support on Eversense DMS Pro, refer to your DMS Pro User Guide, contact your Eversense Representative, or contact Eversense Customer Service at:
1-844-SENSE4U (736-7348)
or support@eversensediabetes.com

Senseonics

© 2026 Senseonics, Inc.
All rights reserved.

Manufactured by Senseonics, Inc.
20451 Seneca Meadows Parkway
Germantown, MD 20876-7005 USA
eversensecgm.com

The Eversense® 365 Continuous Glucose Monitoring (CGM) System is indicated for continually measuring glucose levels for up to one year in people (18 years and older) with diabetes. The system is indicated for use to replace fingerstick blood glucose (BG) measurements for diabetes treatment decisions. Fingerstick BG measurements are required for calibration one time a week after day 13, and when symptoms do not match CGM information or when taking medications of the tetracycline class. The sensor insertion and removal procedures are performed by a health care provider. The Eversense 365 CGM System is a prescription device; patients should talk to their health care provider to learn more.

For safety information, see bit.ly/eversensesafety

Eversense, Eversense 365 Continuous Glucose Monitoring, and the Eversense logo are trademarks of Senseonics, Incorporated. Other brands and their products are trademarks or registered trademarks of their respective holders and are used solely for informative purposes. No relationship or endorsement should be inferred or implied. MKT-002241 REV 01